

COLLIN MALLET DOUBLE HEADS 20 CM 200 GR (STAINLESS STEEL)

Ref. IM01220



10/01/2025 03:58:32 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

COLLIN MALLET DOUBLE HEADS 20 CM 200 GR (STAINLESS STEEL)



Ref. IM01220

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040020372	HS Code	9018.9090
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	4	Bare Unit / Width (cm)	3.4
Packed unit / Height (cm)	2	Bare Unit / Height (cm)	2.8
Packed unit / Length (cm)	20	Bare Unit / Length (cm)	20.7
Packed unit / Weight	0.026	Bare Unit / Weight (kg)	0.026
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

10/01/2025 03:58:32 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr