

## LERICHE FORCEPS CURVED 1X2 TEETH 15 CM

Ref. IP17115



27/08/2024 03:00:35 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France  
[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31  
[info@gsh-med.fr](mailto:info@gsh-med.fr)

# LERICHE FORCEPS CURVED 1X2 TEETH 15 CM



Ref. IP17115

Allows you to perform certain acts in surgery (holding suture threads, compresses, etc.)

## Description

Holtex instruments are all made of top-quality stainless steel (AISI 304, 410, 420 and 420A, depending on the type of instrument) to guarantee good durability and resistance to corrosion.


The steels come into contact with human tissues and must therefore comply with international standards and regulations in force.


Holtex guarantees that its medical devices comply with the 2017/745 EU regulation and ISO 13485.

## Features

CE Marking	CE	Medical device	Class I
Bar code	3661040023120	HS Code	9018.9090
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	15
Packed unit / Weight	0.037	Bare Unit / Weight (kg)	0.037
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

 27/08/2024 03:00:35 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France  
 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31  
 [info@gsh-med.fr](mailto:info@gsh-med.fr)