

LERICHE FORCEPS STRAIGHT 1X2 TEETH 15 CM

Ref. IP17015



 27/08/2024 02:59:55 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

LERICHE FORCEPS STRAIGHT 1X2 TEETH 15 CM



Ref. IP17015

Allows you to perform certain acts in surgery (holding suture threads, compresses, etc.)

Description

Holtex instruments are all made of top-quality stainless steel (AISI 304, 410, 420 and 420A, depending on the type of instrument) to guarantee good durability and resistance to corrosion.

The steels come into contact with human tissues and must therefore comply with international standards and regulations in force.

Holtex guarantees that its medical devices comply with the 2017/745 EU regulation and ISO 13485.

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040023106	HS Code	9018.9090
Warranty (months)	24	Country	Pakistan

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	15
Packed unit / Weight	0.036	Bare Unit / Weight (kg)	0.03

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---

27/08/2024 02:59:55 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr