

WIENER NASAL RAP 22 CM

Ref. IR02422



27/08/2024 05:34:02 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr

WIENER NASAL RAP 22 CM



Ref. IR02422

Features

CE Marking	CE	Medical device	Class Ir
Bar code	3661040029009	HS Code	9018.9090
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	6	Packed unit / Height (cm)	2
Packed unit / Length (cm)	30	Packed unit / Weight	0.062
Bare Unit / Weight (kg)	0.062		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

27/08/2024 05:34:02 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr