

HÉGAR SINGLE ENDED DILATOR STAINLESS STEEL 19 CM N°13



Ref. IS01613



Exemple de bougies Hégar simple, inox.
30 tailles disponibles, nous consulter



27/08/2024 06:04:48 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

HÉGAR SINGLE ENDED DILATOR STAINLESS STEEL 19 CM N°13



Ref. IS01613

Description

Hégar Single Ended Dilator Stainless Steel

30 Sizes available, contact us

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040029832	HS Code	9018.9090
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	1.5
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	0.7
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	19.2
Packed unit / Weight	0.057	Bare Unit / Weight (kg)	0.057
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

 27/08/2024 06:04:48 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

 +33(0)4 42 90 31 31
 info@gsh-med.fr