

HÉGAR SINGLE ENDED DILATOR STAINLESS STEEL 19 CM SET OF 30 N°1 TO 30



Ref. IS01632



27/08/2024 06:15:12 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

+33(0)4 42 90 31 31
 info@gsh-med.fr

HÉGAR SINGLE ENDED DILATOR STAINLESS STEEL 19 CM SET OF 30 N°1 TO 30



Ref. IS01632

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040030029	HS Code	9018.9090
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	2.55
Bare Unit / Weight (kg)	2.4		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

 27/08/2024 06:15:12 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

 +33(0)4 42 90 31 31
 info@gsh-med.fr