

WERTHEIM CLAMP FORCEPS 250 MM

Ref. IP25425



27/08/2024 03:20:37 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

+33(0)4 42 90 31 31
 info@gsh-med.fr

WERTHEIM CLAMP FORCEPS 250 MM



Ref. IP25425

Features

CE Marking	CE	Medical device	Class Ir
Bar code	3661040086217	HS Code	9018.9090
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	8	Bare Unit / Width (cm)	8.6
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	0.2
Packed unit / Length (cm)	26	Bare Unit / Length (cm)	26
Packed unit / Weight	0.106	Bare Unit / Weight (kg)	0.106
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

27/08/2024 03:20:37 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr