

BIPOLAR FORCEPS ANGLED

Ref. IP99998



27/08/2024 05:32:40 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

BIPOLAR FORCEPS ANGLED



Ref. IP99998

Features

CE Marking	CE	Medical device	Class I
Bar code	8964000450154	HS Code	9018908500
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	2	Bare Unit / Width (cm)	2.5
Packed unit / Height (cm)	2	Bare Unit / Height (cm)	2.5
Packed unit / Length (cm)	24	Bare Unit / Length (cm)	24.5
Packed unit / Weight	1.352	Bare Unit / Weight (kg)	1.352
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

27/08/2024 05:32:40 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr