

MOBILE STAND ON WHEELS WITH BASKET FOR DDM TOURNIQUET

Ref. A10701



11/03/2025 23:50:05 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr

MOBILE STAND ON WHEELS WITH BASKET FOR DDM TOURNIQUET

Ref. A10701

The mobile stand for DDM pneumatic tourniquets!

Description

This mobile stand is the perfect accessory to be used with DDM pneumatic tourniquets.

Its mounting plate is offered with a handle facilitating the movements. The pneumatic tourniquet is fixed on the plate with screw-in fixing knobs.

The stand is made of 5 high quality castors from which 2 are equipped with brake. The castors can be locked.

A basket is also available to store cuffs and extension cords.

Made of aluminium and epoxy paint for the gallows and the plate and of powder coated steel for the basket and the base with 5 branches, it is very easy to clean using spray disinfection solution.

Weight : about 15kg

High: 98cm

Features

CE Marking	NONE	Medical device	No
Bar code	3700468201458	Warranty (months)	24
Country	France		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-10	Max temperature	40
Minimum humidity	30	Maximum humidity	40

OTHER FEATURES

11/03/2025 23:50:05 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr