

LUXASCOPE OPHTHALMOSCOPE LED 2.5 V - LUXAMED

LUXAMED®

Ref. F1.416.954



13/03/2025 23:16:50 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

LUXASCOPE OPHTHALMOSCOPE LED 2.5 V - LUXAMED

LUXAMED®

Ref. F1.416.954

Direct ophthalmoscope in pocket size

Description

Direct ophthalmoscope in pocket size

Refraction balance: + 20 bis - 20 D ($\pm$ 1 2 3 4 6 8 10 15 20 D)

5 different apertures: reticle, large circle, small circle, semicircle, green filter

Illumination: LED

Can be dimmed for an individual glare-free brightness

Automatic switchoff after 3 minutes

Life expectancy of LEDs: approx. 50.000 hours by constant current drive

Colour temperature: approx. 4.200 K

Optimized design of head

Two-component handle (aluminium/fibre-glass reinforced), also combinable with LuxaScope LED otoscope- and dermatoscope head

Electronic ON/OFF push-button switch

Power supply: 2x AAA Alkaline batteries, inclusive

Available in :

- White : C1.516.914
- Blue : C1.516.934
- Pink : C1.516.974
- Yellow : C1.516.944
- Kids : C1.916.914
- Black : C1.516.114
- Green : C1.516.954
- Purple : C1.516.964

Features

CE Marking	CE	Medical device	Class I
Bar code	04260381964407	HS Code	90185090
Warranty (months)	24	Country	Germany

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc



SPECIFIC STORAGE CONDITIONS

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr