

## SUCTION UNIT CA-MI ASKIR 36 BR POLYCARBONATE JAR 1L



Ref. BCASPI24



13/03/2025 23:00:19 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France  
[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31  
[info@gsh-med.fr](mailto:info@gsh-med.fr)

## SUCTION UNIT CA-MI ASKIR 36 BR POLYCARBONATE JAR 1L



Ref. BCASPI24

Designed for professional use and equipped with a powerful suction system that can drain up to 36 l/min max, Askir 36BR is lightweight and easy to carry.

### Description

Suction Unit Ca-mi Askir 36 BR, polycarbonate jar 1000 ml

, 36 l / min

Designed for professional use and equipped with a powerful suction system that can drain up to 36 l/min max, Askir 36BR is lightweight and easy to carry.

Main applications: Ambulance car, Emergency Department, Home care, patient transfer, Room, Tracheostomy, Veterinary Med.

- Three different operating options: The AC/DC adapter, rechargeable battery and 12V unit make this surgical Suction Unit a versatile unit guaranteeing CONTINUOUS interventions without the risk of overheating.

- Equipped with a visual and audible alarm system indicating if the battery level is low, the main system is also equipped with a pressure gauge (bar & kPa) and a vacuum regulator.

-

Foot with wheels and optional carrying case

-

Accessories provided: Makrolon 2858 autoclave polycarbonate fluid collection jar, equipped with an anti-reflux valve with different options / 8 x 14 Silicone hoses (autoclavable) and conical connectors o 10-11-12 mm / Antibacterial & Hydrophobic filter (single patient) / Cannula CH20-

Power cord with Schuko plug, Automotive Adapter 12V,

Universal adapter AC / DC-

Motor: Oil-free piston pump requiring no maintenance

- Power supply: 14V 4A with Universal AC / DC adapter / rechargeable Internal battery (12V 4A) Autonomy 60min

- Charging time 240min / Automotive adapter (12V 4A

)-

Maximum vacuum (adjustable): 0.80 Bar -80 kPa -600 mmHg

-

Max free air flow: 36 l/min (without jar)

- Noise level: 63 dBA

 13/03/2025 23:00:19 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31

 [info@gsh-med.fr](mailto:info@gsh-med.fr)

## SUCTION UNIT CA-MI ASKIR 36 BR POLYCARBONATE JAR 1L



Ref. BCASPI24

### Features

CE Marking	CE0123	Medical device	Class IIa
Bar code	8054619100424	HS Code	90189084
Warranty (months)	24	Country	Italy
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	5.5
Bare Unit / Weight (kg)	4.39		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-20	Max temperature	45
Minimum humidity	NA	Maximum humidity	85
OTHER FEATURES			

13/03/2025 23:00:19 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France  
[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31  
[info@gsh-med.fr](mailto:info@gsh-med.fr)