


CUFF AND LOWER LIMB CUFF WITH 2 BLADDERS FOR DDM TOURNIQUET

Ref. **GBS202**

DDM



 18/03/2025 19:09:33 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

 +33(0)4 42 90 31 31
 info@gsh-med.fr

CUFF AND LOWER LIMB CUFF WITH 2 BLADDERS FOR DDM TOURNIQUET



Ref. GBS202

Cuffs and lower limb cuff with 2 bladders for DDM pneumatic tourniquets !

Description

Cuffs and lower limb cuff for DDM pneumatic tourniquets are made of polyester. They are disinfectable and autoclavable and composed of **2 removable silicon bladders** for easy autoclaving.

Those ones are designed in straight form

Several sizes are available

Dimensions of the bladder in cm

Cuff child	2 x (35.5 x 6 cm)
Cuff adult	2 x (59 x 6 cm)
Cuff adult	2 x (69 x 9 cm)

Features

CE Marking	CE	Medical device	Class I
Bar code	3700468203513	HS Code	90189084
Warranty (months)	12	Country	France
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

18/03/2025 19:09:33 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr