

ELITE BAGS - PRACTI - BLUE

Ref. MEPRACT02_B

 **ELITE BAGS**



 12/03/2025 00:03:35 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

ELITE BAGS - PRACTI - BLUE



Ref. MEPRACT02_B

Health assistance bag for home consultations or emergency.

Description

Its excellent internal distribution makes it possible to immediately identify materials.

2 extractable compartments with transparent front part

Isothermal amp maker

External dividers

Inner and outer pockets- Measures: 40 x

30 x 16 cm

Capacity: 19.20 L

Weight: 1.82 kg

Material: 420D polyester

Washable

Maximum recommended capacity: 4 kg

- Sold empty

Features

CE Marking	NONE	Medical device	No
Bar code	8436564542243	HS Code	4202.12.99.90
Warranty (months)	24	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	30
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	16
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	40
Packed unit / Weight	2.05	Bare Unit / Weight (kg)	1.82

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Colors	Blue
--------	------

12/03/2025 00:03:35 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand

13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr