

SPENGLER ARMBAND, VAQUEZ- LAUBRY NANO POCKET



Ref. 330200



12/03/2025 02:10:12 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

+33(0)4 42 90 31 31
 info@gsh-med.fr

SPENGLER ARMBAND, VAQUEZ-LAUBRY NANO POCKET



Ref. 330200

Description

Made of black Oxford-grade nylon, with inner bladder and dual tubing (10 cm & 60 cm), with split leather manometer holder, and Velcro fastening.

Child/Small adult (S), nylon, velcro, black : 330220

Adult (M), nylon, velcro, black : 330210

Adult (L), nylon, velcro, black : 330200

Features

CE Marking	NONE	Medical device	Class I
Bar code	3700446028558	HS Code	9018901000
Warranty (months)	6	Country	Indonesia
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.23
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

12/03/2025 02:10:12 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr