


## LED FIBER OPTIC OTOSCOPE - SMARTLED 5500



Ref. 570517B



17/03/2025 21:05:47 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France  
[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31  
[info@gsh-med.fr](mailto:info@gsh-med.fr)

## LED FIBER OPTIC OTOSCOPE - SMARTLED 5500



Ref. 570517B

The SMARTLED 5500 otoscope has been developed with a lens to allow exceptional diagnosis.

### Description

LED lighting with Optical fibreThe Spengler otoscope SMARTLED 5500® is using LED lighting with optical fibre allowing a direct, homogeneous illumination and without reflections of the auricular canal and the eardrum. Luminance 85K LUX. Bimaterial head in ABS and Stainless steel, screwed head for a safe fixation. Ergonomic ABS grip. End clamp with built-in switch button. Optical fibre for optimal transmission of the light beam. x 3 magnification. Supplied in a case with two LR6 AA batteries, 8 speculums (two of 2.5 mm and four of 4 mm).

Available in 3 colors: black (570517), coral (570517C), blueberry (570517B) and also in a paediatric version (light green 570517P)

### Features

CE Marking	CE	Medical device	Class I
Bar code	3700446041717	HS Code	90189042
Warranty (months)	36	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	12	Packed unit / Height (cm)	4
Packed unit / Length (cm)	20	Packed unit / Weight	0.385
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-20	Max temperature	70
Minimum humidity	10	Maximum humidity	90
OTHER FEATURES			
Colors	Blueberry	Packed unit / Quantity	1

17/03/2025 21:05:47 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France  
[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31  
[info@gsh-med.fr](mailto:info@gsh-med.fr)