

## MALE TUNIC, BOUGAINVILLEA LYOCELL , S0, HOLTEX



Ref. WMBOU07\_95



18/03/2025 19:19:49 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)

## MALE TUNIC, BOUGAINVILLEA LYOCELL , S0, HOLTEX



Ref. **WMBOU07\_95**

### The eco-responsible men's medical tunic

### Description

- White men's tunic with short sleeves,
- Eco-responsible: 50% Lyocell (eucalyptus fiber) and 50% Polyester,
- Officer collar,
- Twill, 215 g / m2,
- 3 plated pockets,
- Buttoning by pressure.
- Available in :  
Size 0 (WMBOU07\_90)  
Size 1 (WMBOU07\_91)  
Size 2 (WMBOU07\_92)  
Size 3 (WMBOU07\_93)  
Size 4 (WMBOU07\_94)  
Size 5 (WMBOU07\_95)

The Holtex eco-friendly and responsible medical tunic brings extra comfort and strength to a classic polycotton tunic: superior quality, increased strength, easy maintenance and maximum ease. Its cut is traditional and modern at the same time. Its officer collar, its short sleeves mounted lapels and its 3 plated pockets are nicely finished with a white border which makes this men's tunic a very neat garment.

Lyocell is a natural material, made from eucalyptus wood pulp.

- Soft, lightweight, comfortable and slightly elastic,
- Easy to maintain, wrinkle resistant, it does not lint.
- Antibacterial, hypoallergenic, breathable, it absorbs moisture (up to 50% more than cotton): it is ideal for sensitive skin.
- Environmentally friendly: Eucalyptus cultivation requires much less water and pesticides than cotton.

### Features

CE Marking	NONE	Medical device	No
Bar code	3661041009529	HS Code	6211433100
Country	Tunisia		

#### DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc

#### SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

#### OTHER FEATURES

Bleaching	Bleaching prohibited	Washing	Wash at 60°C
Model	Male model	Cleaning	Dry cleaning
Packed unit / Quantity	1	Size	Size 5

10/03/2025 15:13:45 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)