

## 3-CHANNEL ECG PAPER, ROLLED

Ref. 260024



 18/03/2025 19:02:22 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France  
 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31  
 [info@gsh-med.fr](mailto:info@gsh-med.fr)

## 3-CHANNEL ECG PAPER, ROLLED



Ref. 260024

### Description

PAPER ROLL dimensions 80MM x 20M FOR CARDIOMATE ECG 3

### Features

|                             |               |                           |            |
|-----------------------------|---------------|---------------------------|------------|
| CE Marking                  | NONE          | Medical device            | No         |
| Bar code                    | 3700446026677 | HS Code                   | 4823400000 |
| Country                     | China         |                           |            |
| DIMENSIONS AND WEIGHT       |               |                           |            |
| Packed unit / Width (cm)    | 5             | Packed unit / Height (cm) | 5          |
| Packed unit / Length (cm)   | 8             | Packed unit / Weight      | 0.102      |
| SPECIFIC STORAGE CONDITIONS |               |                           |            |
| Minimum temperature         | NA            | Max temperature           | NA         |
| Minimum humidity            | NA            | Maximum humidity          | NA         |
| OTHER FEATURES              |               |                           |            |
| Packed unit / Quantity      | 1             |                           |            |

18/03/2025 19:02:22 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France  
[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31  
[info@gsh-med.fr](mailto:info@gsh-med.fr)