

=====

Ref. 338020



 08/03/2025 00:50:13 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

 +33(0)4 42 90 31 31
 info@gsh-med.fr



Ref. 338020

Description

Made of yellow Polymex®, with built-in obturating system to allow decontamination by immersion, and Velcro fastening. Legal markings engraved and non-removable.

Features

CE Marking	NONE	Medical device	Class I
Bar code	3700446028589	HS Code	9018901000
Warranty (months)	6	Country	France
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	9	Packed unit / Height (cm)	1
Packed unit / Length (cm)	13	Packed unit / Weight	0.048
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

08/03/2025 00:50:13 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr