

STETHOSCOPE MAGISTER II DUAL, SILVER

Ref. 507251



10/02/2025 22:54:41 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr

STETHOSCOPE MAGISTER II DUAL, SILVER



Ref. 507251

=====

Description

=====

Features

CE Marking	CE	Medical device	Class I
Bar code	3700446029968	HS Code	90181990
Warranty (months)	60		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	16	Packed unit / Height (cm)	4
Packed unit / Length (cm)	32	Packed unit / Weight	0.43
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-10	Max temperature	40
Minimum humidity	NA	Maximum humidity	95
OTHER FEATURES			
Colors	Gray	Packed unit / Quantity	1

10/02/2025 22:54:41 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr