

CHROME-PLATED LAUBRY STETHOSCOPE / BLACK

Ref. 507300



18/03/2025 19:04:16 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr

CHROME-PLATED LAUBRY STETHOSCOPE / BLACK



Ref. 507300

=====

Description

The famous cardiophone. Single head chestpiece (107 g) made of heavy brass.
Elegant and streamlined for easy handling.
High resolution, fiberglass, acoustic diaphragm, angled binaural and screw-on flexible tips.
Ideal for monitoring blood pressure and general auscultation.

Availbe in black (507300)

Features

CE Marking	CE	Medical device	Class I
Bar code	3700446009151	HS Code	9018.90.20.00-9
Warranty (months)	60		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	16	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	32	Packed unit / Weight	0.393
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Colors	Black	Packed unit / Quantity	1

📅 18/03/2025 19:04:16 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

 +33(0)4 42 90 31 31
 info@gsh-med.fr