

STETHOSCOPE PULSE II SINGLE, BLUEBERRY

Ref. 507700



05/03/2025 20:51:41 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr

STETHOSCOPE PULSE II SINGLE, BLUEBERRY



Ref. 507700

The essential stethoscope. Minimum weight, maximum performance.

Description

HIGH ACOUSTIC QUALITY

- Single tube with high-performance conduit
- High resolution binaural
- High resolution diaphragm

COMFORT

- Anti-cold ring
- Silicone ear tips
- New headphone design more adapted to the anatomy
- Aluminium earhook to reduce weight

Available in Carbon (507702), Blueberry (507700), Coral (507701) or Powder Pink (507703).

Available with double bell (adult bell and bell).

Features

CE Marking	CE	Medical device	Class I
Bar code	3700446029869	HS Code	90181990
Warranty (months)	60		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.226
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Colors	Blueberry	Packed unit / Quantity	1

05/03/2025 20:51:41 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr