

## MOBI WITH ADULT ARMBAND (M) SPENGLER, 4 COLORS TO CHOOSE FROM



Ref. 513700


CE  
0459

DISPOSITIF MEDICAL  
CLASSE  
**Im**  
MEDICAL DEVICE

GARANTIE  
**3**  
3 YEARS  
GUARANTEE



 05/02/2025 22:17:42 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France  
 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31  
 [info@gsh-med.fr](mailto:info@gsh-med.fr)

## MOBI WITH ADULT ARMBAND (M) SPENGLER, 4 COLORS TO CHOOSE FROM



Ref. 513700

**Elegant and colorful, the MOBI blood pressure monitor revolutionizes the aneroid blood pressure monitor thanks to its new patented connectors.**

### Description

The "LIFE LINK®" tubing allows the cuff to be changed in less than a second thanks to its magnetic connectors that ensure the connection and sealing between the pressure gauge and the tubing.

- AMBIDEXTRE : Mobi is specially designed to be ambidextrous, choose your configuration and change the position of the decompression valve, the "LIFE LINK®" connectors and the spoon.

- ERGONOMIC PEAR: The ergonomic pear is suitable for all hand morphologies and offers a volume gain of 11% on traditional pears.

- PRECISION VALVE: MOBI is equipped with a decompression valve whose shape is studied to offer the best sensitivity during decompression but also a very good ergonomics and grip.

- SHOCK RESISTANT: Due to its ABS plastic construction, its double inner shell and its elastomer rear part, MOBI is specially designed to offer great resistance to shocks and falls in everyday life.

- EASY ARMBAND with the revolutionary "LIFE LINK®" connectors, designed in TPU coating to be lighter and easy to clean.

- COLORED: Mobi exists in 4 colors: Carbon, Blueberry, Coral and Powder Pink-

Comes with an M cuff and a LIFE LINK® tubing

- Optional: 5 sizes of armbands (XS, S, M, L, XL)

### Features

CE Marking	CE0459	Medical device	Class Im (measuring function)
Bar code	3700446027957	HS Code	90189010
Warranty (months)	36	Country	France
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.52
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	-20	Max temperature	70
Minimum humidity	20	Maximum humidity	85
<b>OTHER FEATURES</b>			
Colors	Blueberry		

05/02/2025 22:17:42 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France  
[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31  
[info@gsh-med.fr](mailto:info@gsh-med.fr)