

## MOBI MANOPOIRE WITHOUT CUFF SPENGLER



Ref. 513702M



CE  
0459



18/03/2025 19:16:37 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)

# MOBI MANOPOIRE WITHOUT CUFF SPENGLER



Ref. 513702M

**Elegant and colorful, the MOBI blood pressure monitor revolutionizes the aneroid blood pressure monitor thanks to its new patented connectors.**

## Description

Manopaire alone, delivered without cuff and without tubing.

- AMBIDEXTRE : Mobi is specially designed to be ambidextrous, choose your configuration and change the position of the decompression valve, the "LIFE LINK®" connectors and the spoon.

- ERGONOMIC PEAR: The ergonomic pear is suitable for all hand morphologies and offers a volume gain of 11% on traditional pears.

- PRECISION VALVE: MOBI is equipped with a decompression valve whose shape is studied to offer the best sensitivity during decompression but also a very good ergonomics and grip.

- SHOCK RESISTANT: Due to its ABS plastic construction, its double inner shell and its elastomer rear part, MOBI is specially designed to offer great resistance to shocks and falls in everyday life.


- COLORED: Mobi exists in 4 colors: Carbon (513702M), Blueberry (513700M), Coral (513701M) and Powder Pink (513703M)


Note that the Mobi manopaire is to be equipped with a "LIFE LINK®" tubing which allows the cuff change in less than a second and ensures the seal between the manometer and the tubing. This manopaire is to be connected to the easy cuff with "LIFE LINK®" connectors; in TPU coating, it is lighter and very easy to clean.

## Features

CE Marking	CE0459	Medical device	Class Im (measuring function)
Bar code	3700446028268	HS Code	9018901000
Warranty (months)	12	Country	France
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.223
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	-20	Max temperature	70
Minimum humidity	20	Maximum humidity	85
<b>OTHER FEATURES</b>			
Colors	Carbon	Packed unit / Quantity	1

 18/03/2025 19:16:37 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France  
 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31  
 [info@gsh-med.fr](mailto:info@gsh-med.fr)