

MOBI MANOPOIRE WITHOUT CUFF SPENGLER



Ref. 513703M



CE
0459



18/03/2025 19:16:00 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr

MOBI MANOPOIRE WITHOUT CUFF SPENGLER



Ref. 513703M

Elegant and colorful, the MOBI blood pressure monitor revolutionizes the aneroid blood pressure monitor thanks to its new patented connectors.

Description

Manopire alone, delivered without cuff and without tubing.

- AMBIDEXTRE : Mobi is specially designed to be ambidextrous, choose your configuration and change the position of the decompression valve, the "LIFE LINK®" connectors and the spoon.

- ERGONOMIC PEAR: The ergonomic pear is suitable for all hand morphologies and offers a volume gain of 11% on traditional pears.

- PRECISION VALVE: MOBI is equipped with a decompression valve whose shape is studied to offer the best sensitivity during decompression but also a very good ergonomics and grip.

- SHOCK RESISTANT: Due to its ABS plastic construction, its double inner shell and its elastomer rear part, MOBI is specially designed to offer great resistance to shocks and falls in everyday life.

- COLORED: Mobi exists in 4 colors: Carbon (513702M), Blueberry (513700M), Coral (513701M) and Powder Pink (513703M)

Note that the Mobi manopire is to be equipped with a "LIFE LINK®" tubing which allows the cuff change in less than a second and ensures the seal between the manometer and the tubing. This manopire is to be connected to the easy cuff with "LIFE LINK®" connectors; in TPU coating, it is lighter and very easy to clean.

Features

CE Marking	CE0459	Medical device	Class Im (measuring function)
Bar code	3700446028275	HS Code	90189010
Warranty (months)	12	Country	France
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.223
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-20	Max temperature	70
Minimum humidity	20	Maximum humidity	85
OTHER FEATURES			
Colors	Powder pink	Packed unit / Quantity	1

18/03/2025 19:16:00 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr