

VAQUEZ-LAUBRY CLINIC / IN POLYBAG / LARGE ADULT (L)

Ref. 514202



03/03/2025 21:37:24 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

VAQUEZ-LAUBRY CLINIC / IN POLYBAG / LARGE ADULT (L)



Ref. 514202

=====

Description

Cuff-mounted designed to effectively prevent nosocomial infections. The manometer attached on the cuff can be easily removed for liquid decontamination of the cuff in most disinfection solutions on the market. Large, shock-resistant ABS dial and large inflation bulb. 3 cuff sizes available. Comes in a nylon pouch.

Features

CE Marking	CE0459	Medical device	Class Im (measuring function)
Bar code	3700446024758	HS Code	90189010
Warranty (months)	24	Country	France
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.41
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-20	Max temperature	70
Minimum humidity	20	Maximum humidity	85
OTHER FEATURES			
Packed unit / Quantity	1		

03/03/2025 21:37:24 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr