


GREY SPENGLER (S) CUFF (16-24CM) FOR TENSORIC & AUTOTENSIO



Ref. 527462



19/02/2025 21:19:36 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

GREY SPENGLER (S) CUFF (16-24CM) FOR TENSIONIC & AUTOTENSIO



Ref. 527462

Description

Accessory for ES-60® AUTOTENSIO® TENSIONIC® & Effitens®

Features

CE Marking	CE	Medical device	Class I
Bar code	3700446022969	HS Code	90262040
Warranty (months)	6	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.102
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

19/02/2025 21:19:36 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr