


CLAR HEAD MIRROR-ADJUSTABLE HEADBAND FOR ALL SIZES+POUCH



Ref. 570150



19/02/2025 21:08:59 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr

CLAR HEAD MIRROR-ADJUSTABLE HEADBAND FOR ALL SIZES+POUCH



Ref. 570150


Description



Suitable for all head sizes thanks to a flexible metal headband.
Indirect lighting by reflecting mirror.
Suited to lighting required in daily practice.
Built-in battery compartment, two CR 123 batteries
Supplied with carrying bag

Features

CE Marking	CE	Medical device	Class I
Bar code	3700446011390	HS Code	90189084
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.647
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		


19/02/2025 21:08:59 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

 +33(0)4 42 90 31 31
 info@gsh-med.fr