

KERN MGD SCALE, GLASS TRAY, LCD SCREEN, MAX REACH 250 KG

KERN®

Ref. EKMGD02



 12/02/2025 22:30:34 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

KERN MGD SCALE, GLASS TRAY, LCD SCREEN, MAX REACH 250 KG



Ref. EKMGD02

Features

CE Marking	NONE	Medical device	No
Bar code	3361040000199	HS Code	84231090
Warranty (months)	24	Country	Germany
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc
Bare Unit / Weight (kg)	7.5		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			

12/02/2025 22:30:34 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr