


## LED EXAMINATION LAMP WITH 5 FOOT STAND - LUXAMED ®

Ref. U1.161.014



13/03/2025 23:08:07 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)

## LED EXAMINATION LAMP WITH 5 FOOT STAND - LUXAMED ®

Ref. U1.161.014

Optimum lighting for accurate examinations.

### Description

The LED examination light (930mm) is designed to provide optimum illumination using energy-efficient LED technology. It features an aluminium illumination head for efficient thermal convection and heat radiation, with a light intensity of approximately 50,000 Lux at a working distance of 30 cm. The lamp uses a 15 V / 11.25 W LED lamp with 4 high-power LEDs, guaranteeing a life of around 50,000 hours.

The colour temperature is approximately 4,200 K. The gooseneck construction with a 90 cm arm allows the desired position to be maintained.

The lamp has a maximum height of 165 cm and is mounted on a stable metal base with castors, powder-coated to prevent tipping.

It has a multifunctional plug (EU/UK/US).

Available in 6 different colors:

White: U1.111.014

White/ice blue: U1.131.014

White/rapeseed yellow: U1.141.014

White/jade green: U1.151.014

White/delicate lilac: U1.161.014

White/magenta : U1.131.014

### Features

CE Marking	CE	Medical device	Class I
Bar code	04260381961574	HS Code	94052190
Warranty (months)	24		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	5	Max temperature	40
Minimum humidity	10	Maximum humidity	90

#### OTHER FEATURES

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)