


VAQUEZ LAUBRY NANO / VELCRO CUFF / BLACK / CHILD (S)



Ref. 512200



CE
0459



12/03/2025 01:36:18 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr

VAQUEZ LAUBRY NANO / VELCRO CUFF / BLACK / CHILD (S)



Ref. 512200

=====

Description

Lightweight and robust Cuff-mounted, it meets the strictest requirements. Large, lightweight, ABS shock-resistant dial, in accordance with standard ISO 81060-1. Non-deformable Oxford nylon cuff with split leather manometer holder. Large inflation bulb. Supplied in a nylon pouch.

Features

CE Marking	CE0459	Medical device	Class Im (measuring function)
Bar code	3700446007263	HS Code	90189010
Warranty (months)	24	Country	France
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	14	Packed unit / Height (cm)	8
Packed unit / Length (cm)	20	Packed unit / Weight	0.8
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-20	Max temperature	70
Minimum humidity	20	Maximum humidity	85
OTHER FEATURES			
Packed unit / Quantity	1		

12/03/2025 01:36:18 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr