

PNEUMATIC MANUAL TOURNIQUET - SINGLE USE - DDM

Ref. GMBU04



DDM

CE

FABRICATION
FRANÇAISE

DISPOSITIF MEDICAL
CLASSE
I
MEDICAL DEVICE



27/03/2025 21:32:09 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

PNEUMATIC MANUAL TOURNIQUET - SINGLE USE - DDM



Ref. **GMBU04**

Dedicated to emergency situations, the single-use DDM manual pneumatic tourniquet is completely stand-alone, making it a true French innovation!

Description

Designed to accompany the patient from the accident site to the emergency room, DDM's single-use pneumatic tourniquet is accessible and easy to use, thanks to its safety valve.

This single-use pneumatic tourniquet is self-sufficient in its composition: bulb, cuff and valve. The injured patient leaves and keeps the tourniquet on until surgery.

With a wide cuff and foam reinforcement inside, pressure is distributed more widely, avoiding skin lesions and intense pain.

Without a manometer, thanks to a mechanical valve, this manual pneumatic tourniquet avoids overpressure, so pressure does not exceed 300mmHg.

This tourniquet also features a label directly on the cuff for recording patient data, such as patient name, date and time of tourniquet application.

- Available in :

Adult 23-33cm (GMBU04)

Child (GMBU02)

XL Adult 38-50cm (GMCU05)

Features

CE Marking	CE	Medical device	Class I
Bar code	3700468200024	Country	France
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			

📅 27/03/2025 21:32:09 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand

13290 Aix-en-Provence | France

🌐 www.gsh-med.fr

☎ +33(0)4 42 90 31 31

✉ info@gsh-med.fr