

Ref. 332029



28/03/2025 21:13:01 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

+33(0)4 42 90 31 31
 info@gsh-med.fr



Ref. 332029

Description

Cotton cuff with highly resistant rubber bladder.
Twin 7,5 cm tubing fitted with quick CPC coupling. Legal markings printed.
Washable at 30-40°C.

Features

CE Marking	CE	Medical device	Class I
Bar code	3700446008796	HS Code	9018901000
Warranty (months)	6	Country	Indonesia
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.094
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

28/03/2025 21:13:01 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr