

OSCAR 2 ABP MONITOR - SPHYGMOCOR INSIDE

SUNTECH

Ref. DSU230221



CE
0413

DISPOSITIF MEDICAL
CLASSE
IIa
MEDICAL DEVICE

GARANTIE
2 ANS
2 YEARS
GUARANTEE

28/03/2025 21:08:17 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

OSCAR 2 ABP MONITOR - SPHYGMOCOR INSIDE

SUNTECH

Ref. DSU230221

The reference for ambulatory blood pressure monitoring

Description

- Clinically tested by the British Hypertension Society and the European Society of Hypertension, the Oscar 2 M250 offers very high precision measurements- SphygmoCor technology that provides important information on the clinical impact of arterial stiffness and pulse wave reflection, allowing an informed examination of cardiovascular risk
- Extraordinary comfort thanks to the ORBIT cuff
- Silent
- Lightweight and compact
- A mobile application so that the patient can inform his activity via Bluetooth (optional)
- Up to 48 hours of recording
- Analysis and editing software Accuwin Pro 4.0

Features

CE Marking	CE0413	Medical device	Class IIa
Bar code	3661040217536	HS Code	90181990
Warranty (months)	24	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	19	Packed unit / Height (cm)	19
Packed unit / Length (cm)	31	Packed unit / Weight	1.696
Bare Unit / Weight (kg)	0.284		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

28/03/2025 21:08:17 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr