

Ref. 335111



 31/03/2025 21:47:55 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

 +33(0)4 42 90 31 31
 info@gsh-med.fr



Ref. 335111

Description

Cuffs with dual tube bladder & coupling (60 cm tube with Classic® bulb and valve, and 10 cm with male screw-on metal fitting)

Available in :

Adult (M), nylon, velcro, black (335111)

Adult (L), nylon, black (335101)

Features

CE Marking	CE	Medical device	Class I
Bar code	3700446028435	HS Code	9018901000
Warranty (months)	6	Country	Indonesia
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.245
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

31/03/2025 21:47:55 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr