

MAXI +3 WHEEL STAND ADULT CUFF (M)

Ref. 522705



CE
0123

DISPOSITIF MEDICAL
CLASSE
I_m
MEDICAL DEVICE

GARANTIE
2
2 YEARS
GUARANTEE

02/04/2025 22:41:15 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

MAXI +3 WHEEL STAND ADULT CUFF (M)



Ref. 522705

SPENGLER technology in a blood pressure monitor with a dial with a diameter of 147mm.

Description

HIGH MEASUREMENT
ACCURACY• Double tubing for linear pressure rise

COMFORT• Ambidextrous•
Remotely readable measurement•
2.5-meter spiral tubing• Large
volume ergonomic pear for quick inflation•
Storage basket behind the dial• Comes
with cuff wheeled base

• Comes with armband size M, nylon

OPTION
• 3 sizes of cuffs (S, M, L) available

Features

CE Marking	CE0123	Medical device	Class Im (measuring function)
Bar code	3700446028152	HS Code	90189010
Warranty (months)	24	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	4.16
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

02/04/2025 22:41:15 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr