

ARM BLOOD PRESSURE MONITOR SPENGLER TENSIONIC



Ref. 527500



CE
0459



02/04/2025 22:55:52 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr

ARM BLOOD PRESSURE MONITOR SPENGLER TENSIONIC



Ref. 527500

SPENGLER precision in electronic version

Description

The TENSIONIC arm blood pressure monitor is designed for regular therapeutic monitoring. -

-VERY HIGH ACCURACY OF
MEASUREMENTS•

SPENGLER measurement• irregular pulse rate detector•

WHO tricolor risk indicator•

3 average calculations available to monitor the average of morning values

EASY TO

USE• Fully automatic

• Large backlit screen with tricolor WHO risk indicator.

• New armband 22-42 cm size M/L

FAMILY

USE• 2 x 60 time-stamped memories

ADDITIONAL

INFORMATION• Comes with 4 LR6 AA batteries and a carrying case•

Optional AC adapter

• 2 cuff sizes (S, M/L) available as an option

Features

CE Marking	CE0459	Medical device	Class IIa
Bar code	3700446022624	HS Code	90189010
Warranty (months)	36	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.765
Bare Unit / Weight (kg)	0.765		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-20	Max temperature	55
Minimum humidity	15	Maximum humidity	90
OTHER FEATURES			
Packed unit / Quantity	1		

02/04/2025 22:55:52 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand

13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr