

## CA-MI NEW ASKIR 230/12V SUCTION UNIT BR POLYCARB JAR 1L



Ref. BCASPI18



02/04/2025 22:58:38 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)

## CA-MI NEW ASKIR 230/12V SUCTION UNIT BR POLYCARB JAR 1L



Ref. BCASPI18

Electric surgical Suction Unit for draining body fluids and for performing oral, nasal and tracheal aspiration in adults and children. It is ideal for emergency applications and for postoperative treatments at home.

### Description

Three different operating options: equipped with an AC/DC adapter, a rechargeable battery or a voltage of 12V, this surgical Suction Unit is a versatile unit guaranteeing CONTINUOUS interventions without the risk of overheating.

Main applications: Emergency, Home Care, Patient Transfer, Room, Tracheostomy, Med. Veterinary

Equipped with a

- visual and audible alarm system indicating if the battery level is low, the main system is also equipped with a pressure gauge (bar & kPa) and a vacuum regulator.
- Optional roller foot
- Accessories provided: Fluid collection jar with anti-reflux valve with different options / Silicone hoses (autoclavable) 6x10mm and conical connectors o 8-9-10 mm / Antibacterial & hydrophobic filter (single patient) / CH20 cannula / Power cord with Schuko plug / 12V cigarette lighter cable
- Engine: Maintenance-free oil-free
- piston pump3 types of power supply: 14V 4A with Universal Adapter mod. UE60-140429SPA3 (100-240V~ - 50/60Hz - 100VA) / Internal rechargeable battery Pb (12V 4A) - Autonomy 80min - Charging time 240min / Car adapter (12V 1.9A) Maximum vacuum (adjustable): 0.75 Bar - 75 kPa -563 mmHg Max free air flow: 16 l/min Noise level: 63 dBA Operating cycle: Uninterrupted Weight: 4.9 kg /
- 
- 
- 
- 
- Dimensions: 35 x 21 x 18 cm

02/04/2025 22:58:38 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)

## CA-MI NEW ASKIR 230/12V SUCTION UNIT BR POLYCARB JAR 1L



Ref. BCASPI18

### Features

CE Marking	CE0123	Medical device	Class IIa
Bar code	8054619100370	HS Code	90189084
Warranty (months)	24	Country	Italy
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	21	Packed unit / Height (cm)	24
Packed unit / Length (cm)	37	Packed unit / Weight	4.6
Bare Unit / Weight (kg)	3.5		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-25	Max temperature	70
Minimum humidity	1	Maximum humidity	90
OTHER FEATURES			

02/04/2025 22:58:38 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France  
[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31  
[info@gsh-med.fr](mailto:info@gsh-med.fr)