


Ref. C1.426.914




Rechargeable USB



02/04/2025 22:49:58 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr



Ref. C1.426.914

Innovative dermatoscope in pocket size, rechargeable USB

WITHOUT scale.

Description

- Illumination: LED-Ring with 6 LEDs in concentric structure for a highly detailed image, homogenous light
- Maintenance-free, because no lamp replacement necessary
- Removable contact plate for easy cleaning
- Stepless zoom area with 10 x magnification for a distortion-free image
- Optimized design of head
- Can be dimmed
- Automatic switchoff after 3 minutes
- Life expectancy of LEDs: approx. 100.000 hours by constant current drive
- Intensity of illumination: 10.000 Lux
- Colour temperature: approx. 4.200 K
- Two-component handle (aluminium / fibre-glass reinforced), also combinable with LuxaScope LED Otoscope- and Ophthalmoscope head
- Electronic ON/OFF push-button-switch
- Power supply: 1 lithium-ion rechargeable battery 3.7 V, inclusive
- Incl. recharger with USB plug (adapter EU/UK/US)
- Available in white (C1.426.914) and black (C1.426.114)

Features

| | | | |
|------------|----------------|----------------|----------|
| CE Marking | NONE | Medical device | No |
| Car code | 04260381963882 | HS Code | 90189084 |

| | | | |
|------------------|-------|------------------|----|
| Minimum humidity | 10 | Maximum humidity | 90 |
| OTHER FEATURES | | | |
| Colors | White | | |