

## BABY TO TODDLER ELECTRONIC SCALE + BAG

Ref. EHANDI3S

 **Holtex**

CE



 02/04/2025 22:41:53 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31

 [info@gsh-med.fr](mailto:info@gsh-med.fr)

# BABY TO TODDLER ELECTRONIC SCALE + BAG



Ref. EHANDI3S

4-in-1 baby and child scales: ideal for paediatricians and general practitioners

## Description

4-in-1: this scale weighs and measures babies, as well as children thanks to its removable tray and its double horizontal and vertical toise•

3 weighing modes: classic mode, milk absorption and mom + baby weight• Large

capacity up to 100 kg: to weigh mom + baby•

Large ergonomic weighing tray area resistant and optimized to weigh infants (600(W) x 326(W) x 91(H) mm)

• 4 touch keys whose function Manual hold (Hold)•

Hygienic, easy to clean, stable with 4 non-slip rubber feet, easy to carry (Carrying bag sold separately) •

Conversion of units possible in kg/lb and cm/ft•

Large LCD display (79 x 28 mm) with blue backlight, easy to read and automatically rotates 180° when the weight is greater than 20 kg for convenient reading

• Power: 3 x 1.5V AAA batteries (not included)

• Auto On-Off function for power saving

## Features

CE Marking	CE	Medical device	No
Bar code	3361040000021	HS Code	84231090
Warranty (months)	12	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	34	Packed unit / Height (cm)	16
Packed unit / Length (cm)	63	Packed unit / Weight	1.5
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

02/04/2025 22:41:53 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand

13290 Aix-en-Provence | France

 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31

 [info@gsh-med.fr](mailto:info@gsh-med.fr)