



MINI BAG ETRETAT, DE BOISSY, WITH BROWN LEATHER POCKET

Ref. MDETRE02

De Boissy
— DEPUIS 1982



02/04/2025 22:35:06 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr

MINI BAG ETRETAT, DE BOISSY, WITH BROWN LEATHER POCKET

De Boissy
— DEPUIS 1982

Ref. MDETRE02

The mini version of the ETRETAT Medical Bag

Description

Eco-responsible, customizable and "super organized" for a trendy look with all the practical and modern options.

All fabrics, exterior and interior, certified Global Recycled Standard (GRS), have been made from recycled material and remain washable, resistant and water-repellent.

Etretat: it is the gray of the cliffs and the brown of the earth of this emblematic and natural place of Normandy.

If De Boissy wanted to pay tribute to the caregivers and be inspired for his bags in the shape of the traditional leather bag of doctors, the look of Etretat bags contrasts well and truly with the usual classicism. De Boissy played the card of elegance, trend and modernity. And the fact remains that they have been designed to be "ultra practical", to be the best ally of health professionals during their medical visits at home or in community.

The mini version of Etretat and Madrague Medical Bags

- The interior main space also includes 6 pockets, a separator and the 2 mobile compartments: 1 transparent and 1 isothermal, supplied with its reusable cold gel pocket.
- The ventral pocket measures 31x15 cm.
- Always outside: the 3 small zipped pockets, plus the one with elastic for the gel bottle, the detachable-customizable grained leather look pouch, and the external USB port with integrated charging cable.
Same handles, shoulder strap, name holder and non-slip bottom as the large format.
- Dimensions W x H x D: 32 x 28 x 20 cm
- Weight: 1.6 kg
- Maximum recommended content: 4 kg / 15 liters
- Sold empty

The De Boissy bag is put back in its dust bag to protect it when stored. Grey with patterns, this accessory is designed from recycled materials with a cord tightening.



02/04/2025 22:35:06 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

MINI BAG ETRETAT, DE BOISSY, WITH BROWN LEATHER POCKET

De Boissy
— DEPUIS 1982

Ref. MDETRE02

Features

CE Marking	NONE	Medical device	No
Bar code	3661040218977	HS Code	4202220000
Warranty (months)	12	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	54	Bare Unit / Width (cm)	20
Packed unit / Height (cm)	42	Bare Unit / Height (cm)	28
Packed unit / Length (cm)	57	Bare Unit / Length (cm)	32
Packed unit / Weight	1.05	Bare Unit / Weight (kg)	1.05
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

02/04/2025 22:35:06 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr