


## EXAMINATION COUCH STAINLESS STEEL (QUICK ASSEMBLING),BLUE LAVANDER



Ref. OSA003\_L



02/04/2025 22:55:00 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France  
[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31  
[info@gsh-med.fr](mailto:info@gsh-med.fr)

## EXAMINATION COUCH STAINLESS STEEL (QUICK ASSEMBLING),BLUE LAVANDER



Ref. OSA003\_L

### Quick assembly system

### Description

- Structure in square stainless steel tube

-

Adjustable headrest with high resistance automatic rack and pinion

- Maximum load: 150 kg

-

Head and seat padded with very high density foam (30 kg / m<sup>3</sup>), whose thickness is 5 cm, covered with washable plasticized fabric, fireproof (fire resistance classification: M1)

-

Removable roller holder (for rolls 50 cm wide) positionable at the head or feet

- Available in black, Gray and Lavender.

Stainless steel or Epoxy: what's the difference?

Epoxy: the economic choice. The steel is simply coated with a resin to delay corrosion.

Stainless steel: it is steel modified to become completely STAINLESS steel, it does not rust.

And did you know?

Today, more than 90% of stainless steel is made from recycled steel, it's good for the environment!

### Features

CE Marking	CE	Medical device	Class I
Bar code	3661040196152	HS Code	9402900000
Warranty (months)	36	Country	China

#### DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	15	Bare Unit / Width (cm)	63
Packed unit / Height (cm)	68	Bare Unit / Height (cm)	70
Packed unit / Length (cm)	190	Bare Unit / Length (cm)	190
Packed unit / Weight	27.5	Bare Unit / Weight (kg)	24.5

#### SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

#### OTHER FEATURES

Packed unit / Quantity	1
------------------------	---

📅 02/04/2025 22:55:00 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand

13290 Aix-en-Provence | France

🌐 [www.gsh-med.fr](http://www.gsh-med.fr)

☎ +33(0)4 42 90 31 31

✉ [info@gsh-med.fr](mailto:info@gsh-med.fr)