

HE.CO STOP COMPRESSION CUSHION BANDAGE

BIO SYNEX

Ref. SPSECU86



02/04/2025 23:03:29 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr

HE.CO STOP COMPRESSION CUSHION BANDAGE

BIOSYNEX

Ref. SPSECU86

Cushion designed to replace manual lifeguard compression

Description

Allows compression on a wound to stop bleeding
It avoids several manipulations and it frees the rescuer for other interventions.

The elastic band is fixed by a system of Velcro straps.
And a square of foam and a "compress pad" allow the absorption of blood.

Comes with a pair of single-use gloves

Sold individually.

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040116945	HS Code	30059100
Country	China		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	16	Bare Unit / Width (cm)	11.2
Packed unit / Height (cm)	4	Bare Unit / Height (cm)	4.3
Packed unit / Length (cm)	19	Bare Unit / Length (cm)	17
Packed unit / Weight	0.1	Bare Unit / Weight (kg)	0.07
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

02/04/2025 23:03:29 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr