

SPHYGMOMANOMETER EASY 3 BLACK, BLADERLESS CUFF WITH VELCRO ADULT



Ref. DS99421_N



08/04/2025 21:53:21 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr

SPHYGMOMANOMETER EASY 3 BLACK, BLADERLESS CUFF WITH VELCRO ADULT



Ref. DS99421_N

Anti-shock housing

Description

- Blood pressure monitor with cuff without pocket, without latex, washable nylon, black with velcro closure
- Graduation from 0 to 300 mmHg
- Shockproof housing
- Tube length: 65 cm
- Black rubber pear with ergonomic holding spoon
- Comes with nylon storage kit
- Cuff available in infant size (13 to 20 cm, ref DS99521_N) / child (18cm to 26 cm, ref DS99421_N) / adult (26 to 35 cm, ref DS99322_N) / strong arm (32cm to 42 cm, ref DS99603_N) / Multiarmband

Features

CE Marking	CE0123	Medical device	Class Im (measuring function)
Bar code	3661040173030	HS Code	9018901000
Warranty (months)	12	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	12	Packed unit / Height (cm)	8
Packed unit / Length (cm)	18	Packed unit / Weight	0.33
Bare Unit / Weight (kg)	0.5		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-20	Max temperature	70
Minimum humidity	10	Maximum humidity	85
OTHER FEATURES			
Packed unit / Quantity	1		

08/04/2025 21:53:21 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

 +33(0)4 42 90 31 31
 info@gsh-med.fr