

Ref. M25000

DDM

CE

DISPOSITIF MEDICAL  
CLASSE  
**I**  
MEDICAL DEVICE

GARANTIE  
**2**  
2 YEARS  
GUARANTEE



09/04/2025 22:10:50 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)



Ref. M25000

## Description

Pressure cuffs are medical devices designed for use in operating theatres, intensive care units and emergency transport vehicles. They are used to manually accelerate the infusion of blood and solutions when gravity is not sufficient to achieve the desired flow rates.

The GREY FUSE pressure cuff is compatible with all IV and blood bags, and comes in 3 sizes: 500cc (ref. M20500) , 1000cc (ref. M21000) and 3000cc (ref. M25000).

Optimised monitoring and reliability

Reusable and easy to clean/disinfect

Visualisation of the IV bag thanks to a transparent fabric made up of a screen and a polyamide mesh.

Pressure range from 0 to 300mmHg in 10mmHg increments

Pressure gauge accuracy  $\pm 3$  mmHg

Red zone on dial at 300mmHg to prevent overpressure

Pressure relief controlled by a simple tap.

All our cuff are made in our Castelculier factory near Agen.

## Features

CE Marking	CE	Medical device	Class I
Bar code	3700468207276	HS Code	981090
Warranty (months)	24	Country	France
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc
Bare Unit / Weight (kg)	520		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-10	Max temperature	40
Minimum humidity	30	Maximum humidity	40

OTHER FEATURES 0:50 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France  
[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31  
[info@gsh-med.fr](mailto:info@gsh-med.fr)