

## LUXASCOPE AURIS SET LED 2.5 V "COLOUR-EDITION", WHITE / JADE- GREEN

**LUXAMED®**

Ref. AF.416.954



CE



08/04/2025 21:51:51 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)

## LUXASCOPE AURIS SET LED 2.5 V "COLOUR-EDITION", WHITE / JADE- GREEN

**LUXAMED®**

**Ref. AF.416.954**

Complete set for precision diagnosis: Auris LED 2.5V "COLOUR-EDITION", Jade-green otoscope + LuxaScope LED 2.5V "COLOUR-EDITION", Jade-green ophthalmoscope

### Description

Otoscope :

Lighting: LED ring with 6 LEDs (concentric),

no fibre optics, the field of vision is not limited.

- Maintenance-free, as no bulb change is required.

- Possibility of varying brightness intensity.

- Automatic switch-off after 3 minutes

- LED life: approx. 100,000 hours at constant current  
current

- Illuminance: approx. 10,000 Lux

- Colour temperature: approx. 4,200 K

- Optimised head design

- Swivelling magnifier, 3x magnification

- Injection port for pneumatic test

- Handle made from two materials: aluminium and fibreglass

fibreglass, compatible with the head of the Lux dermatoscope and

LuxaScope ophthalmoscope head

- ON/OFF push-button switch

- Power supply: 2x AAA alkaline batteries, included

- Includes 20 grey disposable specula, 10 of which are 2.5mm and 10 of which are  
4.0mm specula

Ophthalmoscope :

- Correction values: + 20 bis - 20 D ( $\pm$  1 2 3 4 6 8 10 15 20 D)


- 5 different filters available: large spot, small spot, green filter,

half-moon and reticle  
08/04/2025 21:51:51 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand

13290 Aix-en-Provence | France

 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31

 [info@gsh-med.fr](mailto:info@gsh-med.fr)

- LED life: approx. 50,000 hours at constant current

current

colour temperature: approx. 4,200°K

optimised head design

Handle made from two materials: aluminium and fibreglass

reinforced, compatible with otoscope and dermatoscope heads

LuxaScope LED

ON/OFF push-button switch

Includes 2 AAA alkaline batteries

GSH

LUXAMED®

LUXASCOPE AURIS SET LED 2.5 V

"COLOUR-EDITION", WHITE / JADE-GREEN

3/3

Ref. AF.416.954

Features

CE Marking	CE	Medical device	Class I
Bar code	04260381964476	Warranty (months)	24
Country	Germany		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-10	Max temperature	45
Minimum humidity	10	Maximum humidity	90
OTHER FEATURES			

08/04/2025 21:51:51 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France  
www.gsh-med.fr

+33(0)4 42 90 31 31  
info@gsh-med.fr