

CARDIOMATE 3 - ELECTROCARDIOGRAPH 3 CHANNELS



Ref. 260003



15/04/2025 21:49:38 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

CARDIOMATE 3 - ELECTROCARDIOGRAPH 3 CHANNELS



Ref. 260003

CARDIOMATE® ECGs are highly reliable instruments.

Description

ECG Cardiomete 3® and 6® are robust and reliable ECG.

Featuring a 7" color screen tilt they allow viewing 12 lead simultaneously.

Simple to use, recording and printing ECG are triggered with only one button.

The Cardiomete 3® & 6® are Cardiomete equipped with the latest technology storage and communication (SD, USB card, network) that can record internally to 10 000 ECG, send them to a computer (optional software) or directly to an external printer.

The Cardiomete 3® & Cardiomete 6® are equipped with automatic interpretation and pacemaker detection.

The internal battery life of 3 hours, their lightness and carrying handle make Cardiomete 3® & Cardiomete 6® the perfect companion in the office or traveling.

- Carrying handheld.

- Tilting screen : 7-inch color LCD display (resolution: 800x480 px)

- SD card included :4GB SD card, 10,000 reports of ECG memory in XML / PDF / BMP / JPG / DICOM

- Autonomy :The operating life is 1080 ECG.

- Automatic interpretation and intervals calculation.

- Pacemaker detection.

- Network connection, USB & printer.

- Printing and rapid analysis.

- Acquisitions and display 12 simultaneous tracks.

- Alphanumeric keypad.

Features

CE Marking	CE0123	Medical device	Class IIa
Bar code	3700446025472	HS Code	9018110000
Warranty (months)	24	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	24.4
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	6.5
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	31
Packed unit / Weight	4.68	Bare Unit / Weight (kg)	2.4
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-20	Max temperature	50
Minimum humidity	25	Maximum humidity	85
OTHER FEATURES			
Packed unit / Quantity	1		

15/04/2025 21:49:38 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr