


SPO2 REUSABLE DIGITAL SENSOR FOR ADULT 3 METERS

SUNTECH

Ref. DSU250899

CE
0413



15/04/2025 21:28:26 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

SPO2 REUSABLE DIGITAL SENSOR FOR ADULT 3 METERS

SUNTECH

Ref. **DSU250899**

Description

SpO2 articulated finger sensor for adult cable with a length of 3 meters for Suntech CT 40 station

Features

CE Marking	CE0413	Medical device	Class IIb
Bar code	0108698870322340	HS Code	90181990
Warranty (months)	3	Country	Turkey
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	16	Packed unit / Height (cm)	3
Packed unit / Length (cm)	13	Packed unit / Weight	0.115
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

 15/04/2025 21:28:26 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

 +33(0)4 42 90 31 31
 info@gsh-med.fr