

LIAN CLASSIC / GREY COTTON CUFF / LARGE ADULT (L)



Ref. 519050

CE
0459

DISPOSITIF MEDICAL
CLASSE
I_m
MEDICAL DEVICE

GARANTIE
3 ANS
3
GUARANTEE



23/04/2025 22:40:15 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr

LIAN CLASSIC / GREY COTTON CUFF / LARGE ADULT (L)



Ref. 519050

=====

Description

Double tubing model with quick coupling. Non-deformable, chromium-plated metal casing with foam manometer cover. The chromium-plated brass valve and double tubing enable highly linear pressure increase and decompression. Rapid cuff replacement thanks to the new CPC connectors. Supplied in a nylon pouch.

CUFF

- Cotton CUFF, velcro fastening
- Size XXS (10-15 cm), XS (15-22 cm), S (20-28 cm, ref : 519029), M (26-34 cm, ref 519012), or L (32-42 cm), S-M-L (ref :519050)
- Available in grey, navy, red or patterns
- Supplied with carrying case

Features

CE Marking	CE0459	Medical device	Class Im (measuring function)
Bar code	3700446024789	HS Code	90189010
Warranty (months)	36	Country	France

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	18	Packed unit / Height (cm)	8
Packed unit / Length (cm)	26	Packed unit / Weight	1.15

SPECIFIC STORAGE CONDITIONS

Minimum temperature	-20	Max temperature	70
Minimum humidity	85	Maximum humidity	85

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---

23/04/2025 22:40:15 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

+33(0)4 42 90 31 31
 info@gsh-med.fr