

LIAN METAL / GREY NYLON CUFF / OBESE (XL)



Ref. 521316



24/04/2025 22:49:35 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr

LIAN METAL / GREY NYLON CUFF / OBESE (XL)



Ref. 521316

=====

Description

Large dial, ambidextrous model (simply rotate to change hands).Metal casing with large dial for easy reading from all angles. Soft PVC manometer cover. Ergonomic spoon handle and inflation bulb. Cuffs with screw-on metal quick coupling. Supplied in a nylon pouch.

Available in :

- Adult (M), cotton, velcro, red (521317)
- Adult (M), cotton, velcro, grey (521311)
- Adult large (L), cotton, velcro, grey (521310)
- Child/Small Adult (S), nylon, velcro, grey (521316)
- Multibrassard (S, M, L) velcro (521330)

Features

CE Marking	CE0459	Medical device	Class Im (measuring function)
Bar code	3700446024871	HS Code	90189010
Warranty (months)	36	Country	France
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.537
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-20	Max temperature	70
Minimum humidity	85	Maximum humidity	85
OTHER FEATURES			
Packed unit / Quantity	1		

📅 24/04/2025 22:49:35 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

 +33(0)4 42 90 31 31
 info@gsh-med.fr