

Ref. AF.416.944



24/04/2025 22:55:41 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

Ref. AF.416.944

Complete set for precision diagnosis: Auris LED 2.5V Blue "colour edition" otoscope + LuxaScope LED 2.5V Blue "colour edition" ophthalmoscope

Description

Otoscope :

Lighting: LED ring with 6 LEDs (concentric),

no fibre optics, the field of vision is not limited.

- Maintenance-free, as no bulb change is required.

- Possibility of varying brightness intensity.

- Automatic switch-off after 3 minutes

- LED life: approx. 100,000 hours at constant current
current

- Illuminance: approx. 10,000 Lux

- Colour temperature: approx. 4,200 K

- Optimised head design

- Swivelling magnifier, 3x magnification

- Injection port for pneumatic test

- Handle made from two materials: aluminium and fibreglass

fibreglass, compatible with the head of the Lux dermatoscope and

LuxaScope ophthalmoscope head

- ON/OFF push-button switch

- Power supply: 2x AAA alkaline batteries, included

- Includes 20 grey disposable specula, 10 of which are 2.5mm and 10 of which are
4.0mm specula

Ophthalmoscope :

- Correction values: + 20 bis - 20 D (\pm 1 2 3 4 6 8 10 15 20 D)

- 5 different filters available: large spot, small spot, green filter,
half-moon and reticle

- Illumination: LED

- Brightness can be dimmed to avoid glare

24/04/2025 22:55:41 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand

13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

- Optimised head design
- Handle made from two materials: aluminium and fibreglass
- Reinforced, compatible with otoscope and dermatoscope heads
- LuxaScope LED

- ON/OFF push-button switch
- Includes 2 AAA alkaline batteries



Ref. AF.416.944

Features

CE Marking	CE	Medical device	Class I
Bar code	04260381964469	Warranty (months)	24
Country	Germany		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-10	Max temperature	45
Minimum humidity	10	Maximum humidity	90

OTHER FEATURES

24/04/2025 22:55:41 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr